MYO-FIT 4

...for your health

ELETTROSTIMOLATORE MUSCOLARE A 4 CANALI
FOUR CHANNEL ELECTRONIC MUSCLE STIMULATOR

FUNCTION
2+2

CE 0123

S. D. T. Co. Ltd.
Shenzhen (PRC)

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
MEDEL MYO-FIT 4
FOUR CHANNEL ELECTRONIC MUSCLE STIMULATOR

Instructions manual

EN

INDEX

STANDARDS	36
CONTRAINDICATIONS	36
INTRODUCTION	37
ELECTRICAL MUSCLE STIMULATION (EMS)	38
TRANSCUTANEOUS ELECTRICAL NERVE STIMULATION (TENS)	38
MEDEL MYO-FIT 4 EQUIPMENT AND ACCESSORIES	39
DESCRIPTION OF DEVICE MEDEL MYO-FIT 4	40
SYMBOLS ON THE MEDEL MYO-FIT 4 DISPLAY	41
PRACTICAL RULES FOR GENERAL USE	42
BEFORE USING MEDEL MYO-FIT 4	42
CONNECTING THE CHARGER	43
POSITIONING OF THE ELECTRODES	43
CARE OF ELECTRODES	46
MAINTENANCE OF THE SELF-ADHESIVE ELECTRODES	47
MAINTENANCE OF THE ELECTRODE CABLE INSULATION	47
SCHEDULING OF THE STIMULATION SESSIONS	48
ALTERNATION OF STIMULATION SESSIONS/VOLUNTARY TRAINING	48
PRELIMINARY SETTINGS: LANGUAGE SETTING	49
USE OF THE MEDEL MYO-FIT 4 STIMULATOR	49
PROGRESSION IN THE LEVELS	52
PROGRAMS	53
ELECTRICITY CONSUMPTION AND RECHARGING	58
SUGGESTIONS FOR BATTERY MAINTENANCE	59
MAINTENANCE INSTRUCTIONS	60
MALFUNCTIONS	60
TROUBLESHOOTING	61
TECHNICAL DATA	62
ELECTRIC STIMULATION	63
SUGGESTIONS FOR A HEALTHY LIFE	63
USEFUL TIPS	64
DEVICE SYMBOLS	65
WARRANTY CERTIFICATE	66

35

The device is based directly on medical technology.

To guarantee your safety, the design, manufacturing and distribution of device are in conformity with the requirements of the European Directive 93/42/CEE.

The device is in conformity with the standard for general safety rules for electro-medical devices IEC 60601-1, the standard for electromagnetic compatibility IEC 60601-1-2 and the standard for special safety rules for nerve and muscle stimulators IEC 60601-2-10.

CONTRAINDICATIONS

Throughout this manual "NOTE" may be found. These Notes are helpful information to aid in the particular area of function being described.


WARNING

- People with pacemakers, intra cardiac defibrillators or other active implants, pregnant women, and children must not use the MEDEL MYO-FIT 4.
- Do not stimulate on the front or sides of the neck, since a drop in blood pressure can occur.
- Do not use the MEDEL MYO-FIT 4 if you are connected to high-frequency surgical equipment, since this may cause burn injuries on the skin under the electrodes, as well as problems with the stimulator.
- Do not use the MEDEL MYO-FIT 4 in the direct vicinity of short-wave or microwave equipment (e.g. mobile phones), since this may affect the output power of the stimulator.
- Keep out of reach from children.

Note

- The stimulator should only be used with skin electrodes intended for nerve and muscle stimulation.
- The electrodes should only be placed on healthy skin. Avoid skin irritation by ensuring good contact between electrodes and skin. Hypersensitivity can occur in rare cases. If skin irritation should occur, treatment with the MEDEL MYO-FIT 4 should be temporarily discontinued. If condition persists, seek medical advice.

EN

- Switch off stimulation before removing the electrodes from the skin. If an electrode comes off, shut off the stimulation before picking it up. Getting electrical stimulation through the fingers is unpleasant but not harmful.
- Exercise caution when using electrodes smaller than 7 cm², as the current density in some cases exceeds 2 mA rms/cm² and increases risk for burn injuries.
- Be careful with stimulating in the heart region.
- CEFAR Medical AB does not take any responsibility for any electrode placements other than recommended.
- Muscle soreness may occur after stimulation with MEDEL MYO-FIT 4, but usually disappears within a week.
- Only Cefar accessories should be used with the MEDEL MYO-FIT 4.
- Use the MEDEL MYO-FIT 4 only as stated in the operating instructions.
- Inspect the equipment prior to use.
- Never use more than one stimulator at the time.
- Never open the battery lid during stimulation.
- The stimulator meets the requirements of MDD 93/42/EEC, including the EMC requirements. This, however, is no guarantee that the equipment will not be affected by electromagnetic disturbances. In case it happens, move away the stimulator from the others devices.
- If the stimulator is not used for some time (approximately 3 months), the non-rechargeable batteries should be removed from the stimulator.

INTRODUCTION

Thank you for choosing the MEDEL MYO-FIT 4. Your new stimulator has been designed with the latest technology in the field of electrical muscle stimulation. Using the MEDEL MYO-FIT 4 regularly will give your training an extra boost, regardless if you are a top athlete, or just exercise to keep fit.

In this manual you are shown how to operate your stimulator safely and how to use it in your training. Be sure to read it carefully before you start training to gain maximum benefits from the MEDEL MYO-FIT 4.

More information on our products and electrical muscle stimulation can be found on our web site: www.medel.it.

ELECTRICAL MUSCLE STIMULATION (EMS)

Electrical Muscle Stimulation (EMS) is successfully used in medical rehabilitation, and to support training in sports at all levels.

The main objective with EMS is to create contractions or vibrations in the muscles. Normally muscular activity is controlled by the central and peripheral nervous systems with electrical signals being sent to the muscles. EMS works in a similar way, but instead the muscle is stimulated with external electrical impulses acting through the skin via the nerve to the muscle. In both cases the muscle reacts and is contracted.

The main benefits of using EMS as a complement to your training:

- During the electrical stimulation, all muscle fibres work simultaneously, making it a very efficient form of training compared to traditional muscle-building exercises.
- You can choose the exact muscle group you wish to train, without the fatigue and energy expenditure resulting from traditional muscle building exercises.
- EMS gives you a faster recovery since the blood circulation is increased and toxin elimination improved.
- EMS provides stimulation of the slow as well as the rapid muscle fibres. The rapid fibres are normally difficult to reach with ordinary training.
- Quick return to activities after injuries.

TRANSCUTANEOUS ELECTRICAL NERVE STIMULATION (TENS)

TENS = Transcutaneous Electrical Nerve Stimulation. TENS is an effective method offering Pain Relief. TENS utilises the nervous system's own pain relief mechanisms to ease both acute and chronic pain.

Good results are observed in acute and long-term pain conditions of many kinds, but above all when the pain originates in joints, the skeleton, muscles, skin, viscera, or nervous system. In some cases, TENS on its own may be sufficient, and in other cases it may provide a valuable complement to other forms of treatment.

With the Pain Relief programs, impulses are sent through the nerves to block the pain impulses. Pain Relief is usually most efficient during stimulation, but the effect can also last after the program has finished. It also increases the

EN

circulation.

The Massage programs increases the circulation and decreases the tension in the muscles.

The Pain Relief programs do not always treat the cause of pain. If the pain persists, seek medical advice.

The programs in the CEFAR EASY can be used on any occasion when Pain Relief or Muscle Relaxation is needed. All the programs can be used as often as required.

Every session should last at least a minimum of 30 minutes up to several hours.


MEDEL MYO-FIT 4 EQUIPMENT AND ACCESSORIES

This device includes:

- A. MEDEL MYO-FIT 4 electro stimulator
- B. 4 cables for electrodes
- C. 4 pairs of electrodes
- D. Instructions manual
- E. Battery pack
- F. Battery charger
- G. Carrying Bag

39


This unit is equipped with multi-function keys. It means that in addition to their primary function, these keys can be used for a whole range of other functions.

1. Turn the device on and off
2. ON/OFF. Language selection and on/off switching
3. Confirm the selection made and enter into the next selection
4. Return to the previous selection
5. Stop stimulation temporarily pressing once to pause and press twice to stop
6. Begin to work
7. Increase or decrease the programs
8. Adjust the intensity of each working channel from 1 mA to 120mA
9. Select Channel A or Channel B
10. Select Channel C or Channel D

40


MEDEL MYO-FIT 4 has four channels and eight electrodes, which mean you can stimulate four muscle groups at the same time with plenty of preset programs. The most innovative features is this muscle stimulation devices that gives you the option of running one program on channels A & B and another program on channels C & D. This is termed "2+2" and permits simultaneous stimulation of different body parts with different programs. For example, you can use a recovery program on a thigh muscle after training while stimulating your neck muscles with a muscle relaxing program. This saves you time and gives more flexibility.

Not only MEDEL MYO-FIT 4 can stimulate the motor muscle of the healthy persons to improve their muscle performance, the device could also stimulate the muscles motor nerves degraded by a process of under- or non-utilization to restore their functional muscular capacities, but under the supervision of a competent health professional.

SYMBOLS ON THE MEDEL MYO-FIT 4 DISPLAY


The use principles presented in this section should be considered general rules. For all programs, it is recommended that you read carefully the information and advice on use presented in this manual. Before a treatment proceed as follows:

- a. Give the device a thorough visual check to see if there are any signs of damage;
- b. Remove all liquid containers from the vicinity of the device because it has not been provided with protection against infiltration by liquids;
- c. Make sure that all the accessories required for stimulation (e.g. electrodes, electrode cables, etc.) are within easy reach.

BEFORE USING MEDEL MYO-FIT 4

1. Connect two electrodes to each cable. Make sure you always insert two electrodes to each electrode cable: one electrode to the positive pole of the cable (red connector) and one electrode to the negative pole (black connector).
2. Apply electrodes to the body. For the correct positioning of electrodes see the table on pages 46 and 50.
3. Connect cables to the MEDEL MYO-FIT 4. The electrode cables themselves are connected to the stimulator through sockets on the top of the device.

! The electrical impulses generated by the device are transmitted to the muscles by self-adhesive electrodes. The choice of the size, the connection and the correct positioning of the electrodes are essential factors for ensuring effective and comfortable stimulation. So they require particular care.

CONNECTING THE CHARGER

EN

The device has considerable operating autonomy. Thanks to the battery pack. However, you do need to recharge it from time to time. To do this, insert the plug of the feeder unit positioned at the rear of the device into the respective socket in the device and then plug the feeder into a power mains socket. It can be recharging when the device is power off.

IMPORTANT: The device cannot work while recharging.

Note: For further information on how to do this, refer to "Electricity consumption and recharging" section of this manual.

When Device is connected to the electrical mains the symbol  remains on constantly to indicate power present.

POSITIONING OF ELECTRODES

P The correct positioning of electrodes is of utmost importance to achieve best results. Before positioning the electrodes on the body part to be stimulated be sure the skin is clean and dry to allow correct conductivity.

Be sure the skin of the area to be stimulated is healthy, without abrasions, wounds, irritations or allergies. Cleanse the skin with water and soap. Do not use alcohol. Remove skin lotions and creams. It might be necessary to epilate very hairy body parts. Do not position electrodes at a distance lower than 3 cm or higher than 30 cm.

Pay utmost attention during stimulation on the cardiac region. Do not stimulate the front and side regions of the neck since this might cause a drop in arterial blood pressure. Small electrodes are more suitable to stimulate small muscles, big electrodes for bigger muscles.

In people with sensitive skin stimulation can cause reddening of the skin under the electrodes. Generally the reddening disappears after about 10-20 minutes.

In any case do not repeat stimulation on the same area until reddening has completely disappeared.

Before applying electrodes to the skin or removing them, be sure stimulation is disabled, that is all channels read 0,0. Do not use the stimulator while driving.


You will find below instructions for the correct positioning of electrodes.

EN

For programs that trigger significant contractions, isometric work is recommended, i.e., the extremities of the limbs should be fixed so that there is no significant movement. Working in this way has a dual purpose: on the one hand, to ensure greater safety and, on the other, to limit shortening of the muscle during the contraction, with the risk of sharp cramps under stimulation.


The following table shows the different muscular groups and gives some very useful indications concerning the best stimulation position to adopt and how to start a contraction voluntarily.

! When removing the electrodes, be sure that the device is turned off, remove the electrodes by pulling them in the direction of the hair growth. To protect the electrode pads, don't pull the lead directly, always reveal from corner of the electrode pads.


EN

EN

⚠ Current international standards require that a warning be given concerning the application of to the thorax: increased risk of cardiac fibrillation.

CARE OF ELECTRODES

Store electrodes on the film side into the sealed bag.

Store electrodes at normal ambient temperature.

In case of loss of adhesion capacity wet the adhesive side with some drops of water and let water be absorbed before positioning electrodes on the film. If the electrodes become gluey or will not stay stuck to the skin despite they have been moistened replace them: electrodes wear out over time and have to be replaced. It is advisable to replace them every 20-40 sessions.

MAINTENANCE OF THE SELF-ADHESIVE ELECTRODES

EN

Proceed as follows to ensure ideal conservation of the electrode layer of gel:

- a. After each use, apply a slight quantity of water on the adhesive gel and let it dry in open air for a few seconds before replacing it on the film.
- b. Tightly close the envelope containing the electrodes to prevent dehydration.
- c. Store the electrodes at room temperature.

It is important to remember that:

- * The duration of the working life of the electrodes depends on their maintenance, the cleanliness of the skin to which they are applied, and the type and intensity of current used.
- * The electrodes must be replaced when they no longer adhere to the skin.
- * The electrodes must be placed only on unbroken skin.
- * Whenever signs of allergy reaction or irritation appear, the stimulation treatment must be suspended and a physician must be contacted.
- * The electrodes have been designed for personal use.

MAINTENANCE OF THE ELECTRODE CABLE INSULATION

Prior to each electro-stimulation session, the state of cables must be checked. If signs of wear such as cuts and/or cracks are observed on their surface, replace the cables with others of the same kind of those originally supplied with the device.

For safety reasons, we recommend never using accessories (such as electrodes, sponges, dials and feeders) other than those provided as standard equipment. The frequency of maintenance, operation control and the check on compliance to EN60601-1 Safety Standards for medical devices to be performed with a secure-tester is annual. The useful working life of the device is guaranteed by the Manufacturer only whenever these maintenance intervals have been duly respected.

NOTE: We recommend having these control operations performed only by the specialized personnel authorized by this latter.

SCHEDULING OF THE STIMULATION SESSIONS

The question of the scheduling of the stimulation sessions during the week only arises in situations where at least two training sessions are to be done in the same week.

For example, anyone who does three sessions a week should do them at the rate of one session every two days (one rest day twice and two rest days once every week). Anyone who does six sessions should do six consecutive days of stimulation with one rest day.

Above seven sessions a week, it is advisable to group several sessions together on the same day to leave yourself one or two complete rest days without stimulation. Anyone who does seven sessions a week should do five days of stimulation, leaving one rest day.

ALTERNATION OF STIMULATION SESSIONS / VOLUNTARY TRAINING

The stimulation sessions can be done outside or during voluntary training. When voluntary training and stimulation are done during the same session, it is generally recommended that the voluntary training be done first followed by the stimulation. This means that the voluntary training is not done on muscle fibres which are already tired. This is particularly important for strength and explosive strength training. However, in resistance training, it can be very useful to proceed in the reverse order. Before the voluntary training, by means of stimulation for resistance, a "specific pre-fatigue session" is carried out on the muscle fibres without general and cardio-vascular fatigue. In this way, the voluntary effort on the "prepared" fibres will push the glycolytic metabolism faster and further.


PRELIMINARY SETTINGS: LANGUAGE SETTING

EN

Before using the unit for the first time, you are strongly advised to take careful note of the counter-indications and safety measures detailed at the beginning of this manual (see chapter I: "Warnings").

The device user can use it with the interface menus displayed in English, French, German, Italian, Spanish or Dutch. In order to select the language to use, press or hold  button for 5 seconds to enter the "language" setting mode, and then choose with  /  language and confirm with ENTER select the language. Press and hold  button to quit for the language setting mode. If no button pressed in 2 minutes, the system will quit the language setting mode automatically.


USE OF THE MEDEL MYO-FIT 4 STIMULATOR

Before using the unit for the first time, you are strongly advised to take careful note of the counter-indications and safety measures detailed at the beginning of this manual.

 The choice of a program and, if necessary, the various stimulation parameters is decisive for ensuring the effectiveness of the desired training or treatment! In this regard, it is particularly useful to consult the chapter "Programs and specific applications" of this manual, which describes each of the programs and presents detailed protocols for the most diverse requirements.


1. Switch the device on
2. Selection of the channel
3. Selection of the program
4. Selection of the treatment position
5. Start the treatment
6. Adjusting stimulation intensity
7. Start/pause/stop the treatment

- 1 MEDEL MYO-FIT 4 is switched On through the button on the device side after through the  button.

NOTA BENE

In case of device inactivity for more than 2 minutes the device will switch off automatically. The stimulator selects automatically the last used program. To change the program follow instructions below.

- 2 Press  button to choose the program for channel A/B channel then press  button to choose the program for C/D channel. The LCD highlights the channel A/B or channel C/D if it is selected by user.

MEDEL MYO-FIT 4 gives your option of running one program on channels A and B and another program on channels C and D.

There are totally 39 preset programs under 7 categories: SPORT, FITNESS, AESTHETIC, VASCULAR, MASSAGE, PAIN, and REHABINATION. Programs of each categories are shown on page 57.


- 3 Press keys  or  to select the program for channels A/B and C/D.

The LCD will display each program name. Press  button to confirm your choice and go for treatment position selection.

- 4 Press  or  to choose cure positions. Press  to confirm your choice. Press  again to go to .

You could select the following treatment position. First press.

Whilst in the *START* phase, for the *Sport* programmes (2 - 11) it is possible to select the level (1-5) using the **P+** and **P-** buttons.

PROGRESSION IN THE LEVELS

In general, it is not advisable to go through the different levels quickly with the intention of reaching level 5 as fast as possible. The different levels correspond to progress in training and you must leave the muscles time to adapt, and time for the overcompensation to be put in place.

The most frequent error is to change from level to level as stimulation occurs at increasingly higher stimulation energies. The number of fibers subjected to stimulation depends on the stimulation energy. The nature and amount of work that these fibers do depend on the program and level. The aim is, first of all, to progress through the electrical stimulation energies and then through the levels. The more numerous the muscle fibers you stimulate, the more numerous will be the fibers that are going to progress. But the speed of progress of these fibers and their aptitude for operating at a higher rating depend on the program and level used, the number of sessions per week and the length of these sessions, and also on intrinsic factors particular to each individual.

The simplest and most usual procedure is to raise the level in the program selected when changing to a new stimulation cycle.


It is also possible to go up a level in the course of the same cycle. In this case, it is advisable not to do so before working a minimum of three weeks at the same level.

Do not change level during an auxiliary or maintenance application. Do not change level either during a short intensive or aggressive cycle of three to four weeks. On the other hand, following the classic approach, during a six-week cycle, after three weeks you can change to the level above. Similarly, in an intensive or aggressive cycle of six to eight weeks, you can go up a level after three to four weeks.

- 5 **Press  to start the treatment.** The treatment timer will start to count down and the output intensity is starting from 0 (no output).
- 6 **Press / to choose the intensity for each channel A\B\C\D.** The intensity range is from 0 to 120.

NOTA BENE: Press  button to choose the intensity for channel A and press again to choose the intensity for channel B. **Then press  button to choose the intensity for channel C and again for channel D.**

During stimulation, the user can vary current intensity by using the / keys corresponding to the channel in use. The value of intensity set is

showed on the display.
 Press ▲ button to increase, 1 each step
 Press ▼ button to increase, 1 each step

7 Press  button to pause the treatment. Press  button again the system returns to the program category mode.

 **ATTENTION:** Whenever you feel anything wrong, press ON/OFF button to turn it off.

PROGRAMS

The device has a considerable number of preset programs in which all the parameters are fixed apart from intensity. The various programs vary in relation to the device and are subdivided into sport, fitness, aesthetic, vascular, massage, pain, and rehabilitation.

SPORT			
N° program	Programs	Uses	Placement of electrodes
P39-1	Potentiation	For optimum muscle preparation immediately before competition.	Depending on the muscle being stimulated, follow the recommended placement. See picture 1-23, 31, 34, 43, 44.
P39-2	Endurance	For athletes who wish to improve their performance in long-duration sporting events.	Depending on the muscle being stimulated, follow the recommended placement. See picture 1-23, 31, 34, 43, 44.
P39-3	Resistance	For competitive athletes who wish to increase their capacity to sustain intense and prolonged effort.	Depending on the muscle being stimulated, follow the recommended placement. See picture 1-23, 31, 34, 43, 44.
P39-4	Strength	For competitive athletes practicing a discipline that requires strength and speed.	Depending on the muscle being stimulated, follow the recommended placement. See picture 1-23, 31, 34, 43, 44.

53

SPORT

N° program	Programs	Uses	Placement of electrodes
P39-5	Explosive strength	For athletes practicing a discipline in which explosive strength is an important factor in performance.	Depending on the muscle being stimulated, follow the recommended placement. See picture 1-23, 31, 34, 43, 44.
P39-6	Fartlek	At the beginning of the season to "re-set" the muscles after a rest period and before more intensive and specific training During the season for those who do not have any preferences for a single kind of performance and prefer to exercise their muscles with different work-out methods.	Depending on the muscle being stimulated, follow commended placement. See picture 1-23, 31, 34, 43, 44.
P39-7	Concentric	Combination of voluntary dynamic concentric contractions with weights and electrically induced contractions.	Depending on the muscle being stimulated, follow the recommended placement. See picture 1-23, 31, 34, 43, 44.
P39-8	Eccentric	Combination of voluntary dynamic eccentric contractions with weights and electrically induced contractions.	Depending on the muscle being stimulated, follow the recommended placement. See picture 1-23, 31, 34, 43, 44.
P39-9	Plyometry	For athletes who on account of injury, must interrupt or limit their voluntary plyometric training sessions For athletes who wish to increase the quantity of plyometric muscular stimulations without risk of injury To reproduce plyometric work on muscle groups for which this type of work is difficult to achieve voluntarily (deltoïd, latissimus dorsi, etc.)	Depending on the muscle being stimulated, follow the recommended placement. See picture 1-23, 31, 34, 43, 44.
P39-10	Hypertrophy	For body-building enthusiasts and users who want to increase their muscle mass For body-builders.	Depending on the muscle being stimulated, follow the recommended placement. See picture 1-23, 31, 34, 43, 44.
P39-11	Stretching	For all athletes who wish to maintain or improve their muscular elasticity To be used after all training sessions or during a specific stretching session.	You have to stimulate the muscle which is antagonist of the stretched muscle Some examples: - Stretching of the calf: see picture 4 - Stretching of the hamstrings: see picture 6 - Stretching of the quadriceps: see picture 5 - Stretching of the biceps. See picture 20.

**EN**

SPORT			
N° program	Programs	Uses	Placement of electrodes
P39-12	Active recovery	To improve and accelerate the muscle recovery after an intensive exercise To be used during the first 3 hours after each session of intensive training or after a competition.	electrodes Depending on the muscle being stimulated, follow the recommended placement. See picture 1-23, 25, 27, 28, 31, 34, 43, 44.
P39-13	Regeneration	To be used the day after competition as recovery training or as a complement to this type of training, the intensity of which can then be reduced.	See pictures no. 1-23 and 25, 27, 28, 31, 34, 43, 44.

FITNESS			
N° program	Programs	Uses	Placement of electrodes
P39-14	Muscle Starter	Improvement of the contractile qualities of insufficiently used muscles; Restoration or improvement of cellular exchanges at the level of stimulated muscles.	Depending on which muscle is to be stimulated, follow the recommended placement see pictures 1-23, 31, 34, 43, 44.
P39-15	Muscle building	For those who wish to gain muscle volume and mass.	Depending on which muscle is to be stimulated, follow the recommended placement see pictures 1-23, 31, 34, 43, 44.
P39-16	Aerobic	Increased muscle ability to sustain long-lasting efforts; improved consumption of oxygen of stimulated muscles.	Depending on which muscle is to be stimulated, follow the recommended placement see pictures 1-23, 31, 34, 43, 44.

AESTHETIC			
N° program	Programs	Uses	Placement of electrodes
P39-17	Muscle tone	To be used to begin with in order to tone up and prepare muscles before more intensive firming work.	Depending on which muscle is to be stimulated, follow the recommended placement, see pictures 1-23, 31, 34, 43, 44. Some examples: Abdominal muscles: picture 10; Buttocks: picture 9; Thighs: pictures 6 and 8; Arms: picture 21.
P39-18	Muscle firming	To be used as the main treatment for muscle firming.	Depending on which muscle is to be stimulated, follow the recommended placement, see pictures 1-23, 31, 34, 43, 44. Some examples: Abdominal muscles: picture 10; Buttocks: picture 9; Thighs: pictures 6 and 8; Arms: picture 21.

AESTHETIC			
N° program	Programs	Uses	Placement of electrodes
P39-19	Shaping	To define and sculpt the body when muscles are already firm.	Depending on which muscle is to be stimulated, follow the recommended placement, see pictures 1-23, 31, 34, 43, 44. Some examples: Buttocks: picture 9; Thighs: pictures 6 and 8; Arms: picture 21.

VASCULAR			
N° program	Programs	Uses	Placement of electrodes
P39-20	Capillarization	These programs are used to increase flow of blood to muscles and when used regularly to increase capillary capacity and promote re-absorption of excess liquids. Particularly indicated for those who engage in endurance sports, during preparation for athletic competition and intense physical work-outs.	Depending on which muscle is to be stimulated, follow the recommended placement, see pictures 1-23, 31, 34, 43, 44.
P39-21	Lymphatic drainage	Deep massage of the stimulated area. Activation of return lymphatic circulation.	Depending on which muscle is to be stimulated, follow the recommended placement, see pictures 1-23, 25, 31, 34, 43, 44.
P39-22	Cramp prevention	To reduce the possibility of cramp happened.	Depending on which muscle is to be stimulated, follow the recommended placement, see pictures 1-23, 25, 31, 34, 43, 44.
P39-23	Heavy legs	Deep massage of stimulated area to reduce heavy feeling of body.	Depending on which muscle is to be stimulated, follow the recommended placement, see picture 25.

MASSAGE			
N° program	Programs	Uses	Placement of electrodes
P39-24	Relaxing massage	Decreased muscular tension. Drainage of toxins responsible for exaggerated increase of muscular.	Depending on which muscle is to be stimulated, follow the recommended placement, see pictures 1-23, 31, 34, 43, 44.
P39-25	Reviving massage	Considerable increase of blood flow in the stimulated region. Improvement of tissue oxygenation; Elimination of free radicals.	Depending on which muscle is to be stimulated, follow the recommended placement, see pictures 1-23, 31, 34, 43, 44.
P39-26	Toning massage	Activation of blood circulation; Recovery of muscular contractile properties; Invigorating effect.	Depending on which muscle is to be stimulated, follow the recommended placement, see pictures 1-23, 31, 34, 43, 44.

**EN**

PAIN			
N° program	Programs	Uses	Placement of electrodes
P39-27	Sophisticated tens	It refers to applications of low-intensity and short duration pulses for pain relief.	On the painful area, try to cover as much as possible with electrodes. See for exemple picture 26.
P39-28	Endorphinic	Block the pain through release of endorphins; Increased blood flow.	On the painful area, try to cover as much as possible with electrodes. See for exemple picture 24.
P39-29	Contracture	Reduced muscular tension has relaxing effect.	On the painful area, try to cover as much as possible with electrodes. See for exemple picture 24.
P39-30	Cervical pain	Block the pain through release of endorphins. Increased arterial follow.	Refer to picture 15.
P39-31	Thoracic back pain	Analgesic current specifically adapted to pain in the thoracic region (below the back of neck and above low back).	Refer to picture 13.
P39-32	Low back pain	Analgesic current specifically adapted to persistent low back pain(lumbar region).	Refer to picture 12.
P39-33	Lumbosciatica	Analgesic current specifically adapted to low back pain (lumbar region) combined with pain in buttocks running down the rear of thigh and sometimes into leg and foot.	Refer to picture 32.
P39-34	Lumbago	Analgesic current specifically adapted to sharp and sudden low back pain (lumbar region).	Refer to picture 33.
P39-35	Arthralgia	Analgesic current specifically adapted to persistent pain in a joint, as in osteoarthritis and rheumatism (all joints in the limbs except the elbow: see the epicondylitis programme below).	On the painful joint, try to cover as much as possible with the electrodes. See for exemple picture 26.
P39-36	Epicondylitis	Analgesic current specifically adapted to persistent pain in the elbow.	Refer to picture 36.

57

REHABILITATION			
N° program	Programs	Uses	Placement of electrodes
P39-37	Disuse-atrophy	It is light with low contraction frequency, to be used to stimulate injured muscles during the early rehabilitation phases.	Depending on which muscle is to be stimulated, follow the recommended placement. See pictures 1-23, 31, 34, 43, 44.
P39-38	Muscle growth	Used primarily to recover muscle functionality after accidents.	Depending on which muscle is to be stimulated, follow the recommended placement. See pictures 1-23, 31, 34, 43, 44.
P39-39	Reinforcement	Used to prior to carrying out aerobic and toning up programs.	Depending on which muscle is to be stimulated, follow the recommended placement. See pictures 1-23, 31, 34, 43, 44.

ELECTRICITY CONSUMPTION AND RECHARGING

Never recharge the stimulator with cables and electrodes still connected to it and in contact with the user and never recharge the battery pack with a charger different from the one supplied with the stimulator!

! WARNING: the device cannot work when connected to the mains.

Whenever the device 4 is being powered by its internal battery, the display's symbol bar will show the following BATTERY "  " symbol.

Key:


-  Battery charged indication
-  Battery run-down indication

We recommend recharging the battery as soon as the battery run-down status is indicated. The device is equipped with an automatic switch-off system that warns the operator with a warning light and eventually switches off the device if it is not recharged in time.

RECHARGING

Proceed as follows to recharge the batteries:

- a. Make sure that the device has been switched off.

- b. Make sure that the device is no longer connected to the patient (the output cables and electrodes must be disconnected).
- c. Connect the feeder to the device. After the battery has been recharged, the display shows the image illustrated up. Leave the device under charge until the display switches off. When fully charged the battery provides about 10 hours of operation.

SUGGESTIONS FOR BATTERY MAINTENANCE

The life of a rechargeable battery depends on the number of recharging/run-down cycles it undergoes and how these cycles are performed. We provide the following suggestions for the longer life of the battery:

Whenever the device is not used frequently, we recommend recharging the battery once a month.

For longer battery life, we recommend discharging the battery as much as possible during battery operation.

EN

MAINTENANCE INSTRUCTIONS

Maintenance and cleaning of the device are extremely easy, provided that the instructions below are complied with:

When not in use store the stimulator and its accessories in the original package in a dry place at a temperature of -20 to $+45^{\circ}\text{C}$.

Avoid contact of the stimulator with water. If necessary, dry the stimulator with a soft cloth.

Cables can be better preserved if they remain connected to the stimulator between sessions.

Do not pull cables or connections.

Do not replace flat batteries with batteries of different type.

MALFUNCTIONS

Should any malfunctions occur while using the device. Check whether the parameters are set to the appropriate form of therapy, adjust the control correctly.

TROUBLESHOOTING

EN

FAULT	ORIGIN	SOLUTION
Batteries are inserted, but nothing appears on the LCD display.	This is not a fault.	Please remove the batteries and replace them after approx. 3 seconds.
The display is working normally, but the user cannot feel any impulses.	Please check whether the plugs on the unit and on the electrode pads are completely inserted.	The plugs must be completely inserted.
The impulse stimulation is weakening.	Ensure that the electrodes are in full contact with the skin. Ensure that there are no foreign particles on the electrodes and that the electrodes are not soiled. Skin impurities.	Press electrodes firmly onto the skin. Clean or replace gel pads. Clean skin surface.
The unit stops during use.	Check status of battery pack.	Change battery pack if necessary.
The skin is turning red and feels like it had been pricked with a needle.	Your skin may be reacting to the gel pads.	Please cancel the session and consult your GP about allergic Reactions.
Even though the unit is set on high, you do not feel any or only a weak stimulation.	Is the skin dirty? Is the electrode pad correctly fastened? Are the batteries empty?	Clean the skin with a damp cloth. Correct the position. Replace the batteries.

If there is any other problem, please return the device to your distributor. Do not try to repair a defective device.

EN

TECHNICAL DATA

MEDEL MYO-FIT 4 is an electrical stimulator, with 4 channels, intended for muscular training and pain relief.

NUMBER OF CHANNELS:	4
PHASE DURATION:	100 - 400 μ s
MAXIMUM FREQUENCY:	1 - 120 Hz
POWER SOURCE:	battery pack 6V 1600 mAh
CURRENT FOR CHANNEL:	0 - 120 mA
EXTERNAL DIMENSIONS:	150 x 90 x 35 mm Weight 410 g (with battery)

REFERENCE TO STANDARDS: Electric Safety Standards CEI EN 60601-1.
Safety Standard IEC 60601-2-10.
Special Rules for nerve and muscle stimulators.
Electromagnetic Compatibility according to CEI EN 60601-1-2.
Class IIa Medical Device according to 93/42/EEC "Medical Device" Directive.

ENVIRONMENT FOR STORAGE AND SHIPPING

Temperature:	-20 °C - +45 °C
Air humidity:	75% max
Air pressure:	700 hPa-1060 hPa

ELECTRIC STIMULATION

EN

Electrostimulation was originally used only in rehabilitation therapy with patients recovering from an injuries or a surgery.

Nowadays it is also used by sporty and common people to tone up or wind down muscles, for beauty or in pain therapy.

It is commonly accepted that electrostimulation can be the least tiring and most efficient method to improve one's physical appearance and to perform high level sport activities in a short time.

Actually, electric stimulation is a very efficient complement in your training, but it cannot replace completely physical activity done through exercising.

Electrostimulation has to be regarded as an "integration" to a precise training schedule. It can both improve muscle quality and act as a real prevention and rehabilitation activity.

Electrostimulation is a special way of training applied to muscles by means of electric impulses (of different frequencies according to the selected program) which contract physiologically the muscle strengthening its action and increasing its mass, just like in ordinary training sessions.

As far as aesthetic is concerned, on many occasions electrostimulation has been wrongly claimed as capable of achieving fabulous results.

It is important to stress that so far there are no scientific studies proving the efficiency to lose weight or eliminate cellulite through electrostimulation.

It is certain that if used with specific programs and accompanied by physical activity and a balanced diet, electrostimulation can contribute significantly to improve physical aesthetic appearance.

Electrostimulation alone cannot work miracles, but it is an important complement to achieve body balance and a sense of wellbeing.


SUGGESTIONS FOR A HEALTHY LIFE

MEDEL MYO-FIT 4 achieves higher efficiency if accompanied by healthy diet and regular physical activity.

A simple daily walk for a couple of miles can help to achieve better results.

Hiking, cycling or swimming are useful aerobic exercises that combined with MEDEL MYO-FIT 4 can help to reach faster results.

Forcing yourself to eat only fruit and salad for a healthy and balanced diet is not necessary, and neither starving to lose weight is recommended. A wider

EN

choice of food is available to keep losing weight, such as potatoes, fish, eggs, chicken, etc.

Low-fat food, meals at regular hours, avoiding high-calory snacks between meals are simple and essential rules.

USEFUL TIPS

- Have your main daily meal before 3.00 p.m., as the peak of metabolism is around 12.00 p.m.
- Eat slowly and chew well to help digestion.
- Avoid alcohol (beer, wine and spirits) due to high-calory content.
- Reduce consumption of tea and coffee because they are stimulant drinks.
- Avoid oils and fat sauces.
- Eat as much fruit and vegetables as possible.
- Avoid salt because it causes water retention.
- Boil, grill or roast food, and avoid to fry.

DEVICE SYMBOLS

EN

Manufacturer

Manufacture Date

Authorized representative in the European Community

P

Atmospheric pressure

RH

Relative Humidity

Temperature restriction

Correct disposal of the product (electric and electronic waste)

(Applicable in the countries of the European Union and in the countries using systems for separate collection of waste)

The symbol on the product or its manuals indicates that the product complies with the directive WEEE (waste electrical and electronic equipment) and must not be disposed of among domestic waste.

At the end of its lifecycle, the user shall dispose of the equipment in specific collection areas. Otherwise, the sanctions foreseen by the current laws about waste are applicable. For detailed information concerning available collection systems, refer to the local waste disposal services.

Type BF applied part

CE 0123

Complying with the community Directive
MDD 93/43/ECC


ATTENTION: Read the instructions

Device connect to the electrical mains.

Warranty conditions

- The device is covered by a 2 YARS warranty from date of purchase for any defects in the design or in the materials used.
- The warranty consists in the replacement and/or repair free of charge of originally defective components.
- The warranty does not cover the accessories supplied with the device and those parts subject to normal wear and tear.
- The device must be repaired by authorized technical service centres only.
- The device must be sent to the authorized service centre for repair within 8 days after the defect is noticed.
- Transport costs for the device are at user's charge.
- Any repair out of warranty conditions are at user's charge.
- The warranty does not cover failures resulting from misuse, abuse of the product or in case the damage is not due to manufacturer's liability (accidental fall, rough transportation, etc.).
- The warranty does not involve any compensation for any kind of damages to persons or things, either direct or indirect, occurred when the device is not operating.
- The warranty is valid from the date of purchase certified by the ticket or the invoice which have to be peremptorily attached to the warranty coupon.
- The lack of warranty coupon properly filled in and validated by a certification of purchase voids the warranty.

COUPON TO BE RETURNED IN CASE OF REPAIRS

**THE WARRANTY IS VALID ONLY IF ACCOMPANIED BY
THE TICKET/INVOICE**

Product type: _____

Model: _____

Serial No.: _____

Date of purchase: _____

BUYER DATA:

Name: _____

Address: _____

Phone number: _____

Defect description: _____

Signature for acceptance of the above mentioned warranty conditions

.....